

WHEN IT'S LIFE OR DEATH, EVERY SECOND COUNTS THINK CELOX

PRODUCT BROCHURE

THE CHALLENGE: SAVING LIVES

FROM POST-TRAUMA BLEEDING

Massive bleeding is the leading cause of preventable death, both in civilian life^{2,3} and on the battlefield.4




THE SOLUTION:

CELOX STOPS BLEEDING FAST

AND HELPS SAVE LIVES

Celox trauma products have helped save the lives of soldiers and civilians since 2006.

- + Celox hemostats use a highquality chitosan formulation, scientifically developed as a fastacting, safe and highly effective haemostat
- + Celox works independently of the body's classic coagulation pathway and does not generate heat⁹
- + Celox hemostats can rapidly promote clot formation, through absorption and dehydration forming a gel-like plug
- + Celox granules work even in challenging cases, including hypothermia and anticoagulation medications, e.g. heparin, warfarin^{9,10}



CELOX RAPID

The most recent addition to the Celox product range is Celox Rapid, our fastest-acting, high-density hemostatic gauze.

- + Celox Rapid contains our latest hemostatic technology with patented Chito-R augmented chitosan, giving fast emergency treatment in critical life-threatening situations
- + Celox Rapid is the UK Ministry of Defence's product of choice¹¹
- + Celox Rapid reduces blood loss by 60% by speeding up packing time and reducing compression time to one minute, or until bleeding stops¹²⁻¹⁴
- + Faster treatment time may improve casualty survival and aid treatment of multiple casualties

PRODUCT BENEFITS:



MODERATE-TO-SEVERE BLEEDING INJURIES + CRITICAL SITUATIONS WHERE COMPRESSION
TIME IS LIMITED + EMERGENCY, PRE-HOSPITAL CARE

HOW CELOX RAPID WORKS

- 1 Celox Rapid contains Chito-R,™ a high-quality chitosan derivative, with the addition of a pharmaceutical-grade bioadhesive
- Chito-R is bonded to a high-density gauze which allows faster packing and greater adherence to surrounding tissue compared to another popular brand. This creates a secure gel plug for maintained hemostasis
- Chito-R activated chitosan granules stick to wet tissue in the wound, improving the dressing's tamponade effect and controlling blood loss¹⁴



PRESENTATION: 1.5m long x 7.6cm wide 'Z-fold' gauze



THICKNESS: 0.68mm

ABSORBENCY: >500%

PRODUCT CODE: FG08839011 (z-fold, 5ft)

ORDERING CELOX PRODUCTS

Celox is sold through a network of distributors. To find out who represents Celox in your area, or to become a distributor, please contact **celoxenquiries@medtrade.co.uk** detailing what country you are in, and with details about your intended use.

Further information

To request a Material Safety Data Sheet (MSDS), or for any other information on the Celox product range, please contact **celoxenquiries@medtrade.co.uk**